Auscultation Sound Library

Stethoscope Quiz

How to use Kikuzo

Auscultation portal site guide

The Sound source posted on the site reproduces a real auscultation sound by listening to the "Kikuzo"

For Auscultation skill acquisition

Supervision

Norio Arakawa

(Former President of the International Society for Hypertension, Professor Emeritus of Fukuoka University)

There is no other way to acquire auscultation skills than to listen to repeated sounds with a stethoscope. Recently, you can learn using a stethoscope recording tool such as a CD, but learning without a stethoscope is not enough. Steso Sound Speaker 3S "Elephant to listen" is an excellent auscultation training device that can realistically reproduce body sounds. In addition, the auscultation portal site has many realistic auscultation sounds. If you use these to practice auscultation, you will be able to learn auscultation techniques more advanced and faster.

About Auscultation portal site

With the advancement and diversification of healthcare, the role of healthcare professionals is being reviewed.

Doctors must acquire the latest knowledge and skills for advanced medical care.

Pharmacists are working hard to acquire physical assessment skills, including auscultation, to play a role in team medicine. Nurses are also becoming more specialized and require more knowledge and skills than ever before.

There is no doubt that a stethoscope is an important clinical tool for healthcare professionals, but nowhere in the field of auscultation training has we found any substantial teaching material.

Therefore, we have developed an original auscultation portal site and a dedicated auscultation speaker that listens to the body sounds such as heart, lungs, and intestines.

The auscultation speaker is developed by a major Japanese acoustic manufacturer, and the auscultation sound source is developed with the support of a specialist. We will work to further expand information in the future.

It is a pleasure to use the listening elephant and the auscultation portal site to improve the auscultation skills and clinical knowledge of everyone who is engaged in medical treatment and all those who aim for medical treatment.

Kiyoshi Fujiki, Representative

Director Telemedica Inc.

Auscultation portal site TOP screen

<https://3sportal.telemedica.co.jp>

local=
staging=
session_lifetime=720
★ LoginAt= Now=2020-03-25 20:58:53

EN
Login

Learn auscultation sounds! Portal site


3&S
stetho-speaker Kikuzo

Home

Kikuzo library

Stethoscope quiz

Video

Auscultation is acquired quickly

Kikuzo to a stethoscope speaker

Auscultation master

To be grateful
"Thank you for your consultation!"

The Special Event

Realistic sound auscultation app

Kikuzo Library

Stethoscope quiz

Video physical assessment

Listen repeatedly using a stethoscope

Supervision
Norio Arakawa
Former President of

There is no other way to acquire auscultation skills than to listen to repeated sounds with a stethoscope. Recently, you can learn using a stethoscope recording tool such

News Letter Vol.5

Writing
Mr. Shoji Kudo
(President, Tuberculosis Prevention Association, Professor Emeritus, Nissai Medical)

Basics of Physical Assessment "Lung Auscultation" Lung auscultation, brought on by Laennec in France at the beginning of the 19th century, is the basis of physical

Purchase

1. When listening to the auscultation sound, Auscultation speaker [Kikuzo] Or use earphones
2. Kikuzo [Recruiting facility] Is [here](#).
3. Posted sound source : 411 Cases ([Sound source list](#))

Auscultation member library

Member registration is [here](#).

PR [Introduction]
Register for the auscultation member library and listen to "heart sounds" for free!

PR [Introduction]
Register for the auscultation member library and listen to "lung sounds" for free!

HEART	Stepwise 0 Introduction	more >		00:00 / 00:00
HEART	Stepwise 1 Sound split	more >		00:00 / 00:00
HEART	Stepwise 2 S2 split 1	more >		00:00 / 00:00
HEART	Stepwise 3 S2 split 2	more >		00:00 / 00:00
HEART	Stepwise 4 Atrial sound 1	more >		00:00 / 00:00
HEART	Stepwise 5 Atrial sound 2	more >		00:00 / 00:00
HEART	Stepwise 6 Ejection sound 1	more >		00:00 / 00:00

Received "EdTech Special Category Award" at "Japan e-Learning Grand Prize".

* EdTech is a coined word that fuses Education x Technology (education x technology).

Speaker for auscultation training
[Kikuzo]

Realistic auscultation sounds in the low range that could not be realized until now can be reproduced!

[See details.](#)

[Click here for usage.](#)

Detailed Information Screen :

Information can be checked while listening to Auscultation sound

- ・ 「  Press ">" to play the auscultation sound.
- ・ 「 **more** 」 Press to display detailed information.


Heart sounds

more >

 |  00:28/01:06

IV sound (S4gallop, benper tuning) _002

IV 音 (excessive heart sounds) ・ S4 gallop / hypertrophic cardiomyopathy (HCM), dilated cardiomyopathy (DCM), etc. Caused by obstruction of rapid inflow / contraction of left atrium / audible as "wa-ka-ta" / IV sound occurs immediately before I sound and is masked by I sound making it difficult to hear / squeezing stethoscope Then, you can not hear the sound, so you can use it to distinguish it from the sound split / listen on the bell surface / # cardi_sound4_002_60bpm

Auscultation sound name:	IV sound (S4 gallop, tuned horse tuning) _002
Auscultation sound type:	Abnormal sound
On behalf of the disease:	Abnormal left ventricular myocardium (hypertrophic cardiomyopathy, dilated cardiomyopathy, etc.)
Auscultation site:	Apex

Detailed Information : Slide unfolds in a schematic diagram of auscultation sound

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
00:28/01:06

IVsound (S4gallop tuned horse) _002

IV sound (excessive heart sounds) · S4 gallop / hypertrophic cardiomyopathy (HCM), dilated cardiomyopathy (DCM), etc. Caused by obstruction of rapid inflow / contraction of left atrium / audible as "wa-ka-ta" / IV sound occurs immediately before I sound and is masked by I sound making it difficult to hear / squeezing stethoscope Then, you can not hear the sound, so you can use it to distinguish it from the sound split / listen on the bell surface / # cardi_sound4_002_60bpm


Auscultation sound name:	IV sound (S4 gallop, tuned horse tuning) _002
Auscultation sound type :	Abnormal sound
On behalf of the disease:	Abnormal left ventricular myocardium (hypertrophic cardiomyopathy, dilated cardiomyopathy, etc.)
Auscultation site:	Apex

IV sound (S4gallop, benper tuning) _002

IV sound (excessive heart sounds) ・ S4 gallop / hypertrophic cardiomyopathy (HCM), dilated cardiomyopathy (DCM), etc. Caused by obstruction of rapid inflow / contraction of left atrium / audible as "wa-ka-ta" / IV sound occurs immediately before I sound and is masked by I sound making it difficult to hear / squeezing stethoscope Then, you can not hear the sound, so you can use it to distinguish it from the sound split / listen on the bell surface / # cardi_sound4_002_60bpm

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
Auscultation sound name: IV sound (S4 gallop, tuned horse tuning) _002

Auscultation sound type: Abnormal sound

On behalf of the disease: Abnormal left ventricular myocardium (hypertrophic cardiomyopathy)

Auscultation site: Apex

Heart Sound

You can audition on the auscultation portal site TOP screen

00:00/01:00


Physiological Sound

Physiological sound (excessive heart sound) / Listened to in healthy young people (three years or younger) with thin chest wall and good myocardial extensibility / Sound of blood flowing from the left atrium bumping into the apex (ventricular filling sound) / You can hear "nat-to-ku" ./ There are physiological sounds and pathological sounds in the sound. Scrutinize as a target sound / Listen to the left side and listen on the bell /

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
Auscultation sound name: Physiological sound (ventricular filling sound) _002

Auscultation sound type: Normal sound

On behalf of the disease: Healthy young people

Auscultation site: Apex


Mid-diastolic noise (diastolic rumble, left ventricular filling sound)

Mid-diastolic murmur (diastolic rumble, left ventricular filling sound), mitral valve opening sound (OS), hypertonia / mitral valve stenosis / stenosis delays mitral valve opening, diastolic (Ⅱ音 and Ⅰ音)
Between: Listen to mitral valve opening sound (OS) at mitral valve opening time) / Listen to middle diastolic rumble (left ventricular filling sound) / Listen to bell on left side lying position / # cardi_mitral_stenosis_001

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


Auscultation sound name: Mid-diastolic murmur (diastolic rumble, left ventricular filling sound), mitral valve opening sound _001

Auscultation sound type: Abnormal sound

On behalf of the disease: Mitral stenosis (MS)

Auscultation site: Apex

Systolic ejection noise (organic noise), arterial ejection noise


Systolic ejection noise, arterial ejection noise / aortic valve stenosis / If there is a valve stenosis, listen for ejection noise (produced by rapid ejection of blood to artery) immediately before systolic noise / Since the sound is heard 0.04 to 0.1 seconds later than the sound, it sounds like a sound split / After the ejection sound, the ejection noise (gradually increasing and decreasing type) becomes the sound of the sound (half moon valve closing time) Prolonged to just before / Ejective murmur is caused by pressure difference between ventricle and artery when there is resistance to blood flow near meniscal valve (aortic valve / pulmonary valve) / Subvalvular stenosis (hypertrophic obstruction) Cardiomyopathy) does not produce an ejection sound and can be distinguished from valve stenosis / #

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

Auscultation sound type: Abnormal sound

On behalf of the disease: Meniscal (aortic, pulmonary) stenosis (AS, PS)

Auscultation site: Aortic valve area


Normal heart sounds_apex_002

Normal heart sound / I sound lasts longer than the II sound / The shorter the interval between the next sound is, the main vibration of the I sound / I sound is the atrioventricular valve (mitral valve / tricuspid valve) closing time (QRS wave Of the aortic valve, which constitutes the end part of the sound / The sound is the vibration at which the blood hits the orifice when the semi-lunar valves (aortic and pulmonary valves) are closed / the apex (TM area) 領域 sound> II sound / # cardi_normal_base_002

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


Heart Sound Library

> Heart sounds posted on the auscultation portal site

- Normal heart sound
- I Muted
- I Hypertone
- II Split tone
- II Increased aorta component
- II Hyperpneumo pulmonary component
- Physiological sound (ventricular filling sound)
- Pathological sound (S3 gallop,奔馬音)
- IV Sound (S4 gallop, 調馬音)
- Overlapping horse tune (gallop)
- Four-part tuning
- Mitral valve opening sound
- Mid systolic click
- Systolic ejection noise (functional noise)
- Systolic ejection noise (organic noise), arterial ejection noise
- Musical noise / systolic ejection noise
- Mid-systolic ejection noise
- Systolic reflux noise
- Mid-diastolic noise (diastolic rumble, left ventricular filling sound)
- Diastolic reflux noise
- Ventricular extraperitoneal contraction (PVC, VPC)
- Sinus tachycardia
- Sinus bradycardia
- Continuity noise

[More than 60 sound sources](#)

※The stethoscope library is updated as needed. Please confirm the latest information on auscultation portal site.

Lung Sound

You can audition on the auscultation portal site TOP screen


00:00/01:00

Sound rhonchi_004 / Katan storage

Sputum retention / Inhalation / expiration / Vibration caused by viscous secretions sticking to the airway wall, or due to vibration of the narrowed airway wall due to airflow / Snoring sound "Gugu" / Inhalation Listening to the low-pitched continuous ra (rhonchi) from the first half / Listening to the slightly low-pitched continuous ra (rhonchi) during expiration / # respiratory_rhonchi_004

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
Auscultation sound name: Sound rhonchi_004 / Katan storage

Auscultation sound type: Abnormal sound

On behalf of the disease: Sputum retention, bronchiectasis, bronchial asthma, COPD, chronic bronchitis, bacterial pneumonia, heart failure

Auscultation site: Propagation on chest wall close to sound source / widespread


Sound crackles_005 / interstitial lung disease

Interstitial lung disease / Inspiration / Sound of sudden opening of the peripheral airway that was obstructed during expiration / Large intermittent buzzing sound from the beginning of the inspiration phase / # respiratory_finecrackles_005

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



Auscultation sound name: Sound crackles_005 / interstitial lung disease

Auscultation sound type: Abnormal sound

On behalf of the disease: Interstitial lung disease, early stage of pulmonary edema, non-typed pneumonia

Auscultation site: Dorsal and precordial lower lung fields / widespread

Lung Sound

You can audition on the auscultation
portal site TOP screen

00:00/01:00


Alveolar breath sounds_002

Alveolar breathing sound / Low soft soft sound / Volume is inspiration > expiration / inspiration time:
expiration time = 1: 2, but after mid-expiration the volume is low and cannot be heard / no pause
(pause) / # respiratory_vesicular_002

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
Auscultation sound type: Normal Sound

On behalf of the disease: —

Auscultation site: The back and precordium in contact with the peripheral lung

Lung Sound Library

> Lung sounds posted on the auscultation
portal site

- Tracheal breathing sounds
- Bronchoalveolar breath sounds
- Alveolar breath sounds
- Twisting sound (fine crackles)
- Blisters (coarse crackles)
- Flute sound (wheezes)
- Chirp (rhonchi)
- Squawk (squawk)
- The snowfall sound (snowball)

[And more than 35 sound sources](#)

※The stethoscope library is updated as needed. Please confirm the latest information on auscultation portal site.

Bowel Sound

You can audition on the auscultation portal site TOP screen

00:00/01:00


Intestinal sound

Transverse colon site / Glue reduction / 14 seconds, listening for gurus at around 59 seconds / Peristaltic sound propagates throughout the abdomen, so auscultation points are sufficient at one or two places / Normal: 5 times or more per minute (5 to 15 (Every second) Hearing sound can be heard, Decrease: 1 to 2 times a minute (constipation, peritonitis, paralytic ileus, intestinal dysmotility due to surgery, disappearance: No gleaming sound for 5 minutes (functional ileus, peritonitis), hyperactivity : Always audible (infectious enteritis, diarrhea, ileus)

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
Bowel Sound and Other

> Bowel sounds and other information posted on the auscultation portal site


Other Sound

- Bowel sound
- Mechanical bowel obstruction (metallic sound)
- Thyroid rales
- Vascular murmur
- Dialysis shunt sound

[And more than 20 sound sources](#)

※The stethoscope library is updated as needed. Please confirm the latest information on auscultation portal site.

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

Auscultation Quiz Library

1. Lung sounds Quiz_1 Guess the correct lung sounds (5 questions)
2. Heart sounds Quiz_1 Guess the correct heart sounds (6 questions)
3. Lung Sound Quiz_2 Which sound is suspected of interstitial lung disease? (6 questions)
4. Heart sounds Quiz_2 Heart failure due to aortic stenosis (AS) is increasing with aging. Early detection possible by auscultation
5. Lung sounds Quiz_3 Do you understand these lung sounds? (6 questions)
6. Heart sound Quiz_3 Is this heart sound normal? Abnormal? (5 questions)
7. Lung sounds Quiz_4 Do you understand these lung sounds? (6 questions)
8. Heart sound Quiz_4 Do you understand this heart sound? (7 questions)
9. Lung sound basic Quiz_1 Let's compare two sounds (6 questions)
10. Heart sound basic Quiz_1 Let's assign the listening part of heart sound! (10 questions)

※Auscultation quiz is updated from time to time. Please confirm the latest information on auscultation portal site.

Telemedica Inc.

〒227-0055 9-9 Azumigaoka, Aoba-ku, Yokohama-shi, Kanagawa

Emai 3sp@telemedica.co.jp

URL <https://3sportal.telemedica.co.jp>